

## 0.2 ml RE-Spec (50-100 mesh) Columns for Whole Rock Bulk REE Separations

### Column Procedure:

- |                        |  |
|------------------------|--|
| 1. clean resin         | 3 ml 0.05 N HNO <sub>3</sub>             |
| 2. precondition column | 3 ml 2 N HNO <sub>3</sub>                |
| 3. load sample         | dissolved in 0.2 ml 2 N HNO <sub>3</sub> |
| 4. wash                | 0.1 ml 2 N HNO <sub>3</sub>              |
| 5. wash                | 0.1 ml 2 N HNO <sub>3</sub>              |
| 6. wash                | 0.1 ml 2 N HNO <sub>3</sub>              |
| 7. wash                | 0.5 ml 2 N HNO <sub>3</sub>              |
| 8. wash                | 0.5 ml 2 N HNO <sub>3</sub>              |
| 9. wash                | 1.0 ml 2 N HNO <sub>3</sub>              |
| 10. wash               | 1.5 ml 2 N HNO <sub>3</sub>              |
| 11. elute REE          | 4 ml warm* 0.05 N HNO <sub>3</sub>       |

\* use of warm acid helps overcome a kinetic effect which results in severe tailing of the REE when room temperature acid is used.

### Resin Cleaning:

RE-Spec resin is a rare earth element selective extraction chromatographic material produced by Eichrom Industries, Darien, IL. We preclean the resin by placing an aliquot in a small Teflon bottle which is then filled with quartz-distilled water, shaken well, and placed on a hotplate set at around 50°C overnight. The next day the water is decanted and the procedure is repeated ten times. Once the resin is ready for use, it is transferred to a dropper bottle where it is dispensed as needed. The resin in the columns is disposed of after every use.

For details on this resin, see:

Horwitz, E.P., Chiarizia, R., Dietz, M.L., and Diamond, H., 1993. Separation and preconcentration of actinides from acidic media by extraction chromatography. *Anal. Chim. Acta*, 281: 361-372. (note: this paper discusses TRU-Spec, an analogous extraction chromatographic material to RE-Spec, but with a lower ratio of CMPO/TBP)

### Column Construction and Cleaning:

Start with 5.5 cm KORVEX 4 × ½ shrink tubing. Fit the shrink tubing over the end of a standard test tube with approximately 2.5 cm of the tubing extending beyond the end of the test tube. Shrink this portion of the tubing with a hot air gun as far as it will go. This will result in a column with an interior diameter of around 3.6 mm. Place a polystyrene frit in the end of the column and adjust it so that the volume of the column is 0.2 ml. The reservoir (the unshrunk portion of the column) has a volume of 3 ml.

After use, the columns are placed in a Teflon bottle which is then filled with clean water and placed on a hotplate at a low setting for a couple of days. After this, the water is decanted, the columns are rinsed, and 7 N HNO<sub>3</sub> is added to the bottle and placed back on the hotplate at a low setting for a couple more days. After this, the HNO<sub>3</sub> is decanted and the columns are rinsed with clean water and are ready for use.

Questions? Contact Todd Housh (housh@mail.utexas.edu)